

*QMSA Letters, 3(2025)*

Spin-layers and packing of large spin-networks

Reino Laatikainen , Pekka Laatikainen and Henri Martonen  
Spin Discoveries Ltd., Kuopio, Finland

Last update Jan. 10<sup>th</sup>, 2025

Eight spins, Yess!



LAOCOON group, which  
gave the name for the early  
QMSA program

3 Layers, yeah!

From LAOCOON to  
ChemAdder



Greetings  
from  
Kuopio,  
home of  
ChemAdder

# Spin-layers


## How to simulate the pentyl CH<sub>3</sub> signal

- The A<sub>3</sub>-particle is strongly coupled to BB' particle, which leads to the characteristic triplet. - The BB' particle forms here the 1<sup>st</sup> Spin-Layer.
- The A<sub>3</sub>-particle couples to CC' particle with <sup>4</sup>J couplings which are not usually resolved with normal line-width: - The CC' particle forms here the 2<sup>nd</sup> Spin-Layer.
- The A<sub>3</sub>-particle couplings to DD' particle are very small and may lead at most to small broadening of triplet lines, DD' forms the 3<sup>rd</sup> Spin-Layer... However, DD' is coupled strongly to CC' which is coupled to BB', and if the chemical shift differences are small (as they typically are for (CH<sub>2</sub>)<sub>n</sub>-systems), also the 2<sup>nd</sup> and 3<sup>rd</sup> layer particles affect the outlook of the A<sub>3</sub>-triplet, through *the second-order effects*.

Structure:

System:

Layer:



To get nearly complete description for the CH<sub>3</sub> signal, one needs the 3 layers, and to optimize the BCD chemical shifts! ChemAdder builds up automatically the sub-systems when the number of layers is given. The sub-systems giving the BB', CC', DD' and EE' transitions are built up in the same way .. See the octyl simulation on next page.



# Spin-system packing

- Compression of spin-networks

- $XX'RR'KK'BB'CC'DD'EE'H_3 \rightarrow \underline{XX'RR'KK'}BB'CC' + KK'\underline{BB'CC'DD'EE}H_3 + DD'EE'\underline{H_3}$ . The underlined species are got from the system.
- Transitions separated by  $< 0.01$  Hz and belonging to the same species are combined, if also their derivatives are similar.
- 58 000 000  $\rightarrow$  24 000 Transition, from 64 to 4 sec

- Multithreading

- Analysis of 16 serum samples: 16 min  $\rightarrow$  4 min.



"Spin dust"


*The  $H_3$ -signal is composed of thousands of non-degenerate transitions – which yield its diagnostic outlook.*

*The effect is not rare, for example leucine methyl signal has a similar shape. See more at page 10.*


***How to describe the shape in another way than QMSA, we ask?!***

# Large spin-networks: Cholesterol

*All the protons coupled with each other through the **spin-network***



**H3, H6 and Me18 out of this range:**




46 protons → 36-spin particles → 14 sub-systems (not all shown above) → 368001 transitions → 9322 lines.

Simulation time < 1 sec, with different line-width for each spin-particle.

# D<sub>3</sub>-vitamin 80 MHz simulation<sup>1</sup> quality and time depend on the *number of layers (1-3)*

35





## 36.

n

---



# FELIX\* of MILLIONS TRANSITIONS

**Experimental:** normal  $^1\text{H}$  spectrum,

*Spin-system of 74 protons :*

XX'RR'AA'BB'CC'DD'EE'  $\text{H}_3$  (Ar-octyl) +


XX'RR'ZZ' BB'CC'H $_3$  (O-hexyl) +

XX'RR'ZZ' BB'CC'DD'EE'H $_3$  (O-octyl) +

XX'RR'ZZ'BB'(CC') $_3$ DD'EE'H $_3$  (O-decyl) +

XX'YY'ZZ'

(5 mol% chiral ingredient with F was not included in QMSA)




Non-zero intensities: ca. 1 750 000 (93% of theor. intensity), Lines after packing: ca. 135 000

*Although the  $\text{CH}_2$ -shift order is somewhat unclear, the analysis gives the alkyl lengths with a fair confidence*

\* S.K.Ahola, L.P.Ingman, R.Laatikainen, J.Sinkkonen & J.P.Jokisaari,  $^{21}\text{Ne}$  and  $^{131}\text{Xe}$  NMR study of electric field gradients and multinuclear NMR study of the composition of a ferroelectric liquid crystal. *J.Chem.Phys.* **149**, 234901 (2018); <https://doi.org/10.1063/1.5052499>



## The FELIX methyl signal carries the information from the alkyl populations



The Layers is an easy way to define the level of approximation and optimize the calculation time

## Full MENU

1. Analyse the system at highest approximation level: 1 layer and broadening.
2. Remove the broadening and optimize system with 1-2 layers.
3. Finalize with 3 layers.