

```
&CASEINFO: %P1=human, %P2=male, %age=50, &type=5, &class=12
&DEFAULTS:
    ORIGINAL = C:\URINE\1_URINE220218.JDX
    BATCHFILE = C:\CHEMADDER\SCRIPTS\URINE.BAT
    PROFILE = C:\CHEMADDER\PROFILE\SEARCHPROFILE.TXT
```

```
&INC C:\CHEMADDER\URINE\HIPPURATE.ASL      POPULATION= 0.82 ***
&INC C:\CHEMADDER\URINE\CITRATE.ASL        POPULATION= 0.72 ***
&INC C:\CHEMADDER\URINE\GLYCINE.ASL        POPULATION= 0.38 ***
&INC C:\CHEMADDER\URINE\ME3-AMINE-OXIDE.ASL POPULATION= 0.32 ***
&INC C:\CHEMADDER\URINE\TAURINE.ASL        POPULATION= 0.29 ***
&INC C:\CHEMADDER\URINE\CYSTEINE.ASL       POPULATION= 0.23 ***
&INC C:\CHEMADDER\URINE\CREATINE.ASL       POPULATION= 0.16 ***
&INC C:\CHEMADDER\URINE\HISTIDINE.ASL      POPULATION= 0.15 ***
&INC C:\CHEMADDER\URINE\GLYCOLATE.ASL      POPULATION= 0.15 ***
&INC C:\CHEMADDER\URINE\ISOCITRATE.ASL     POPULATION= 0.17
&INC C:\CHEMADDER\URINE\GLUCOSE.ASL        POPULATION= 0.13 ***
&INC C:\CHEMADDER\URINE\GLUTAMINE.ASL      POPULATION= 0.13 ***
&INC C:\CHEMADDER\URINE\ETHANOLAMINE.ASL   POPULATION= 0.13 ***
&INC C:\CHEMADDER\URINE\ARABINITOL.ASL     POPULATION= 0.11 **
&INC C:\CHEMADDER\URINE\MANNITOL.ASL       POPULATION= 0.11
&INC C:\CHEMADDER\URINE\TRIGONELLINE.ASL   POPULATION= 0.11 ***
...
&INC C:\CHEMADDER\URINE\ALPHA-KETOISOVALERATE.ASL POPULATION= 0.01 *
&INC C:\CHEMADDER\URINE\SYRINGATE.ASL      POPULATION= 0.01
&INC C:\CHEMADDER\URINE\N-ACETYLPUTRESCINE.ASL POPULATION= 0.01
&INC C:\CHEMADDER\URINE\13-DINH2-PROPANE.ASL POPULATION= 0.01
&INC C:\CHEMADDER\URINE\TRANS-FERULATE.ASL POPULATION= 0.01
&INC C:\CHEMADDER\URINE\MALEATE.ASL        POPULATION= 0.0
```

&END of FILE

<http://chemadder.com>

HOLISTIC QMSA - A NEW DIMENSION OF QMSA:

R. Laatikainen and P. Laatikainen

SPIN DISCOVERIES Ltd, Kuopio, Finland

HOLISTIC qQMSA

A few clicks from a set of data to a presentation!

- qQMSA = quantitative Quantum Mechanical Spectral Analysis
- ASL = Adaptive Spectral Library: NMR Spectra are stored as spectral parameters, so that the spectrum (since it has been analyzed once) can be simulated in any field


HOLISTIC qQMSA: combined analysis of N spectra yields more & better information than N independent analyses for the same spectra !

- **The fundamental PROBLEM of qQMSA:** chemical shifts depend on sample !
 - Ranges typically < 0.01 ppm, but even 0.10 ppm for some common metabolites.
- **The SOLUTION:** *PatternSearch* and *PriorKnowledge* algorithms (=FZZY *Search*), are based on the fact that the variations are not independent but they correlate.
- The more spectra have been analysed, the better correlations => **HOLISTICS**
 - Batch analysis, < 2-5 min/sample

A few clicks from spectra to diagram !

See video: <http://www.chemadder.com/software.html>

HOLISTICS: simultaneous analysis of N spectra yields more & better information than N independent analyses for the same spectra



The more spectra analysed, the better model!


```

AVERAGE.PMR - Notepad
File Edit Format View Help
&SpinAdder SEARCH file
*ChemAdder2017.05 C:\CHEMADDER\JDXLIB\AVERAGE.PMR
%SEARCH FILE

&CONTROL PARAMETERS & FIT INFO
  ORIGINAL = C:\CHEMADDER\JDXLIB\AVERAGE.QMT
  SPECTRUM = C:\CHEMADDER\JDXLIB\AVERAGE.QMT
  BATCHFILE = C:\CHEMADDER\JDXLIB\JDXLIB.BAT
  PROFILE = C:\CHEMADDER\PROFILE\SEARCHPROFILE.txt
  TEMPLATE = C:\CHEMADDER\JDXLIB\AVERAGE.PMR
  FIELD = 600.402822
POINT RESOLUTION = 0.137549
SCALING FACTOR = 1.000000
REFERENCE MMOL = 9.0000
CONCENTRATION = 42.1216
  LINE WIDTH = 3.300
  GAUSSIAN = 0.000
  RRMS = 0.143
  R-FACTOR(%) = 98.522

&INC C:\CHEMADDER\ASLIBS\METLIB\tsp.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\lactate.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\acetate.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\butyrate.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\valeric.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\PROPIO.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\etoh.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\GLUCOSE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\sucrose.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\23bud.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\MEOH.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\HCOOH.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\MANNITOL.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\MANNOSE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\XYLOSE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\GALACTOSE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\SORBITOL.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\13pdio1.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\PHENALN.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\TYROSINE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\HISTIDIN.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\SUCCINICACID.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\ACETONE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\ALANINE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\valine.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\ISOLEUCINE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\LEUCINE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\ISOBUTYRATE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\XMETHYL.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\MEMALONICACID.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\ARGININE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\ASPARAGN.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\ASPARTAT.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\BETAINE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\buoh.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\CITRATE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\creatine.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\CRN.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\GLUTAMIN.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\glyce.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\LYSINE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\PROLINE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\SERINE.ASL
&INC C:\CHEMADDER\ASLIBS\METLIB\THREONIN.ASL

  POPULATION= 2.420
  POPULATION= 13.082
  POPULATION= 8.235
  POPULATION= 5.457
  POPULATION= 0.250
  POPULATION= 1.845
  POPULATION= 2.175
  POPULATION= 5.339
  POPULATION= 1.912
  POPULATION= 0.939
  POPULATION= 0.877
  POPULATION= 1.039
  POPULATION= 3.291
  POPULATION= 0.096
  POPULATION= 2.345
  POPULATION= 0.180
  POPULATION= 0.512
  POPULATION= 5.122
  POPULATION= 0.262
  POPULATION= 0.294
  POPULATION= 0.209
  POPULATION= 0.559
  POPULATION= 1.270
  POPULATION= 2.697
  POPULATION= 0.633
  POPULATION= 0.980
  POPULATION= 0.339
  POPULATION= 1.268
  POPULATION= 1.388
  POPULATION= 0.178
  POPULATION= 4.471
  POPULATION= 1.675
  POPULATION= 0.470
  POPULATION= 1.914
  POPULATION= 1.914
  POPULATION= 0.729
  POPULATION= 0.653
  POPULATION= 2.086
  POPULATION= 1.948
  POPULATION= 4.641
  POPULATION= 1.914
  POPULATION= 1.492
  POPULATION= 7.593
  POPULATION= 3.310


&CONSTRAINTS (in equations X0 = 1.0)
  IGNORE(PPM): 5.15 to 4.70

&END OF FILE

```

RECIPE (PMR-file) after FZZY SEARCH

HOLISTIC CYCLE of bioproduct analysis with *ChemAdder*



Under testing ...



<http://chemadder.com>